



Rapid Vascular On-Growth with Titanium Nanotube vs. Acid Etched Surfaces

Authors:

David Detwiler, PhD

Sabrina Huang, PhD

Alan Kraft, BSNE

Rapid Vascular On-Growth with Titanium Nanotube vs. Acid Etched Surfaces

Authors: David Detwiler, PhD; Sabrina Huang, PhD; Alan Kraft, BSNE

Introduction

Nanotube technology has been around for a while and well-studied by multiple labs around the world. Khosravi et al. has done an animal study looking at the vascularization around implants¹. Utilizing a literal glass window, they monitored vascular on-growth and vessel maturation that led to new bone formation around the implant. The test materials were micron-roughened acid-etched titanium (TiMA) and an anodized titanium nanotube surface (TiNT). The authors found that the nanotube surface TiNT significantly outperformed the micron-roughened acid-etched surface TiMA in speed of vascularization, bone volume around implant and bone to implant contact.

The foundation of creating new bone is the vasculature that will feed the new bone, provide oxygen and remove waste^{2,3}. Without the resources needed to make bone, the body will make fibrous tissue that can't stabilize orthopedic, trauma or dental implants. With insufficient resources, the body will heal more slowly.

Titanium dioxide nanotubes have demonstrated the ability to interact with the body to encourage a healing response that further recruits vasculature^{1,2}. Nanotube surfaces with specific size and range of size features on titanium implants enable these effects through increased protein binding and spacing of cell attachments.^{4,5}

Methods

The work was performed by Khosravi et al. Implants composed of grade 4 titanium were machined into a cross-shaped structure, 4 mm in diameter, with a central hole 2 mm in diameter. Flutes were machined to limit the bone contact at implantation and encourage vascular growth and leave a void space that can be filled by bone. The titanium surfaces were shot-peened and acid-etched to create a nanostructure similar to dental implants (TiMA). Half of the implants were further modified by anodization to create titanium dioxide nanotubes with an outer diameter of about 70 nm (TiNT).

Titanium cranial implants were implanted into the skull bone of mice. The schematic of the implantation is shown in Figure 1. A glass window was glued over the implant to prevent infections and allow for visualization of the vascular growth and maturation process. To visualize the vascular system, fluorescent dye was injected into the tail vein of mice prior to visualization of vascular structures with confocal microscopy. The glass window allowed for sequential visualization of the same implant/mouse and the tracking of vascular growth onto the implants. After 42 days, the animals were sacrificed, and the implants were scanned using micro-CT imaging to quantify bone growth onto the implants.




Figure 1: Placement of the titanium implants into the skull of the mouse allows bone contact. Covering the wound with a glass cover allows for visualization of the vascular bed with fluorescent microscopy and tracking of vascular on-growth and maturation. This image is modified from the work of Khosravi et al.




Fig. 2. In vivo longitudinal microvascular response to TiMA and TiNT cranial implants. Representative images demonstrating formation and development of the peri-implant neovascular network from Day 3 to Day 42 around the TiMA surface and the TiNT surfaced implants. Scale bars: 500 μ m. Images are stacks of tiled scans of the entire craniotomy at the maximum intensity projection; the depth of the field is 0.5 mm, which is equal to the thickness of the implant. This image is modified from the work of Khosravi et al.

Results

The vascular growth onto the titanium implants with nanotubes (TiNT) occurred more quickly than on the micron-roughened acid-etched surface (TiMA). At each time point there was more vascular volume on the TiNT surface than the TiMA surface. Vessel coverage was also complete at an earlier timepoint (Day 15) on nanotubes, as shown in Figure 2. Hierarchical vessel maturation was also evident on the nanotube surface at Day 11, with continued maturation out to 42 days. The TiMA surface showed delayed vessel maturation relative to the nanotube surfaces. There was also more bone-to-implant contact 76% vs. 11% and greater bone volume 35% vs. 10% with the nanotube surface vs. acid etching alone, graphs available in full publication¹.

Discussion

Nanotube surfaces have been shown to modify the *in vivo* behavior of multiple tissue types, including vasculature and bone¹⁻⁷. There are multiple studies by various peer reviewed groups that have shown the benefits of nanotube surfaces *in vitro*⁴⁻⁷. But *in vitro* studies can't replicate the complexity of what happens *in vivo*. This study allows a literal window into the process of vascular growth onto implants that results in new bone formation. Improvements in vascularization with nanotube surfaces on implants will directly benefit patients by allowing more complete osseointegration.

Conclusion

The authors of this publication demonstrated that nanotube modified surfaces (TiNT) have a significant advantage over micron-roughened acid-etched surfaces (TiMA) *in vivo*. Nanotubes are better at recruiting vasculature sooner and growing bone more quickly. Nanotube modified titanium surfaces, such as nanoVIS Ti™ Surface Technology, provide solutions for implant manufacturers looking for a permanent surface that can produce better implant integration and outcomes for patients.

References

1. Khosravi, N., Maeda, A., DaCosta, R. S. & Davies, J. E. Nanosurfaces modulate the mechanism of peri-implant endosseous healing by regulating neovascular morphogenesis. *Commun Biol* **1**, 72, doi:10.1038/s42003-018-0074-y (2018).
2. Wu, B., Tang, Y., Wang, K., Zhou, X. & Xiang, L. Nanostructured Titanium Implant Surface Facilitating Osseointegration from Protein Adsorption to Osteogenesis: The Example of TiO₂. *Int J Nanomedicine* **17**, 1865-1879, doi:10.2147/IJN.S362720 (2022).
3. Filipowska, J., Tomaszewski, K. A., Niedźwiedzki, Ł., Walocha, J. A. & Niedźwiedzki, T. The role of vasculature in bone development, regeneration and proper systemic functioning. *Angiogenesis* **20**, 291-302, doi:10.1007/s10456-017-9541-1 (2017).
4. Oh, S. *et al.* Stem cell fate dictated solely by altered nanotube dimension. *Proc Natl Acad Sci U S A* **106**, 2130-2135, doi:10.1073/pnas.0813200106 (2009).
5. Rosa, A. L. *et al.* Nanotopography Drives Stem Cell Fate Toward Osteoblast Differentiation Through $\alpha 1\beta 1$ Integrin Signaling Pathway. *Journal of Cellular Biochemistry* **115**, 540-548, doi:10.1002/jcb.24688 (2014).
6. Chamberlain, L. Macrophage Inflammatory Response to TiO₂ Nanotube Surfaces. **2**, doi:10.4236/jbnb.2011.23036 (2011).
7. Yu, W. P. *et al.* Titanium dioxide nanotubes promote M2 polarization by inhibiting macrophage glycolysis and ultimately accelerate endothelialization. *Immun Inflamm Dis* **9**, 746-757, doi:10.1002/iid3.429 (2021).