



Titanium Nanotube Surfaces Improve Osseointegration in a Diabetic Rat Model

Authors:

David Detwiler, PhD
Sabrina Huang, PhD
Alan Kraft, BSNE
Kreigh Williams, BS



Titanium Nanotube Surfaces Improve Osseointegration in a Diabetic Rat Model

Authors: David Detwiler, PhD; Sabrina Huang, PhD; Alan Kraft, BSNE; Kreigh Williams, BS

Background

Diabetes is highly prevalent in the United States and around the world. In the United States, there are an estimated 38.4 million people (11.6%) with diabetes and 97.6 million people (38%) with pre-diabetes¹. For those people 65 years or older, 48.8% have pre-diabetes¹. Diabetes is a chronic progressive inflammatory disease with effects on multiple systems including vasculature, nervous and skeletal systems^{2,3}. The inflammatory state increases bone resorption while limiting bone deposition, leading to osteoporosis over time^{2,3}. Bone loss can lead to fractures that require orthopedic implants⁴. Impaired healing due to diabetes can lead to complications such as implant loosening, device failure, and infections⁴.

Titanium nanotubes have shown the ability to improve outcomes in osteoporotic models due to hormone dysregulation⁵. In this paper, Zhang et al. explored the effect of titanium nanotubes on osseointegration under diabetic conditions using a rat model, where diabetes-associated inflammatory and impaired vascularization typically inhibit implant integration⁶.

Methods

Machined titanium screws were used as the control group, while the experimental group (TNTs) consisted of machined titanium screws that were further anodized to create the nanotube surface. The diabetic animal model was created by injecting the 9-week-old rats with streptozotocin to kill off the beta-islet cells in the pancreas. The rats' body weight, intraperitoneal glucose tolerance test (IPGTT), and insulin tolerance test (ITT) were evaluated to confirm the development of type 2 diabetes mellitus (T2DM).

Screws with and without nanotubes were placed into the distal femur of the rats. The study also included control rats and T2DM rats without implants. Micro-CT scanning was performed to record bone growth around the implants at 4-, 8-, and 12-week timepoints to generate data on bone volume/total volume (BV/TV), trabecular thickness (Tb.Th), trabecular number (Tb.N), bone surface/bone volume (BS/BV), and trabecular separation (Tb.Sp).

Results and Discussion

Micro-CT imaging showed obvious osteoporosis in the distal femurs of the T2DM rats compared with the control rats, as shown in Figure 1. Consistent with these findings, the T2DM rats resulted in significant bone loss in the distal femur and dramatically reduced weight gain compared with the control rats (weight data are not shown here but are reported in the original publication).

At 4-, 8-, and 12-week timepoints, micro-CT images showed increased bone growth around the nanotube surface relative to the control surface, as shown in Figure 2.

Quantification of the bone structure at 12 weeks showed greater bone volume/total volume (BV/TV), thicker trabecular structures (Tb.Th), and a higher trabecular number (Tb.N) for the nanotube surface than for the control surface, as listed in Table 1. In contrast, bone surface/bone volume (BS/BV) and trabecular separation (Tb.Sp) were lower for the nanotube surface. Due to the increased trabecular number and thickness, the spacing between trabeculae on the nanotube implants was lower, resulting in higher bone volume.




Figure 1: Micro-CT images of distal rat femur with bone loss due to T2DM and normal control at 12-weeks post injection. This image is modified from the work of Zhang et al⁶.




Figure 2: Micro-CT images of bone around control titanium and titanium nanotube surfaces screws at 4-, 8-, and 12-weeks after implantation. This image is modified from the work of Zhang et al⁶.

Table 1: Micro-CT quantification of bone around implants at 12 weeks after implantation. This table is modified from the work of Zhang et al⁶.

Group	BV/TV	Tb.Th	Tb.N	BS/BV	Tb.Sp
Control	21.7 ± 2.6	28.3 ± 4.7	2.1 ± 0.2	9.6 ± 0.4	423.2 ± 23.5
TNTs	25.2 ± 4.1	46.7 ± 6.6**	3.5 ± 0.2**	8.9 ± 1.2	354.1 ± 21.3**

**p < 0.01

Conclusion

Diabetes is a pervasive, chronic and progressive disease state that inherently causes bone loss, leading to vascular, neurological and skeletal complications¹⁻³. Fractures are common as the disease progresses into osteoporosis, requiring orthopedic implants or amputation if the implants fail to achieve proper healing¹⁻³. Titanium nanotubes present a technology upgrade for the orthopedic implants that doesn't rely on drugs or supplements but rather works with the host biology to lower inflammation, recruit vasculature, improve stem cell differentiation to heal and retain bone around the implant⁷⁻⁹. In this study, the combination of higher BV/TV, Tb.Th, and Tb.N, together with lower BS/BV and Tb.Sp, indicated more robust bone formation for the nanotube surface than for the control surface under diabetic conditions in the rat model.

References

1. C. f. D. Control. (2024), vol. 2024, pp. Statistics for diabetes in the United States of America based on 2021 study results.
2. V. Sundararaghavan, M. M. Mazur, B. Evans, J. Liu, N. A. Ebraheim, Diabetes and bone health: latest evidence and clinical implications. *Ther Adv Musculoskelet Dis* **9**, 67-74 (2017).
3. B. Wu et al., A narrative review of diabetic bone disease: Characteristics, pathogenesis, and treatment. *Front Endocrinol (Lausanne)* **13**, 1052592 (2022).
4. A. Arteaga et al., Diabetes as a Risk Factor for Orthopedic Implant Surface Performance: A Retrieval and. *J Bio Tribocorros* **7**, (2021).
5. H. Kang et al., Titania-Nanotube-Coated Titanium Substrates Promote Osteogenesis and Suppress Osteoclastogenesis via Integrin $\alpha v \beta 3$. *ACS Appl Bio Mater* **5**, 5832-5843 (2022).
6. Q. Zhang et al., Improved osseointegration of dental titanium implants by TiO_2 nanotube arrays with self-assembled recombinant IGF-1 in type 2 diabetes mellitus rat model. *Nanotechnology Reviews* **12**, (2023).
7. W. P. Yu et al., Titanium dioxide nanotubes promote M2 polarization by inhibiting macrophage glycolysis and ultimately accelerate endothelialization. *Immun Inflamm Dis* **9**, 746-757 (2021).
8. N. Khosravi, A. Maeda, R. S. DaCosta, J. E. Davies, Nanosurfaces modulate the mechanism of peri-implant endosseous healing by regulating neovascular morphogenesis. *Commun Biol* **1**, 72 (2018).
9. S. Oh et al., Stem cell fate dictated solely by altered nanotube dimension. *Proc Natl Acad Sci U S A* **106**, 2130-2135 (2009).